

## Understanding the Patient Tab

The **Imaging Devices** drop down menu allows you to select different devices, like **EVA** or a **camera**.

These 4 buttons allow you to **create a new patient**, **anonymous patient**, **clear** the selected patient, or **print**

File Edit **Image Devices** Tools Preview Help

**Patient** Acquire Chart View **Diagnosis View**  
Monroe, Marilyn

Current Patient: Marilyn Monroe (No study selected).

Patient Information Results of Last Search



- Patient/Study
- Search Results
  - Chanel, Coco
  - Monroe, Marilyn
    - Study - Nov 9, 1973 7:25 AM
    - Study - Oct 21, 1976 4:22 AM
    - Study - Sep 30, 1964 11:01 PM
  - Pantomime, Cephalo
  - Sinatra, Frank

| Patient/Study | Patient Id      | DOB          | Gender | Last Study       | # Studies | Dentist          |
|---------------|-----------------|--------------|--------|------------------|-----------|------------------|
| 512985698745  | 123456789ABCDEF | Dec 28, 1950 | F      | 3/15/11 12:21 PM | 4         | Einstein, Albert |
| 365236523     | 3654AFDGT       | Jul 29, 1953 | F      | 3/4/11 11:10 AM  | 2         | Einstein, Albert |
| 3654AFDGT     |                 | Dec 26, 1978 | M      | 3/10/11 12:19 PM | 3         | Einstein, Albert |

Search Options:

Quick Search

DICOM Worklist Patients

Deleted Patients Only

Flagged Patients

Anonymous Patients Only

Studies Over Last Month

Yesterday/Today Studies

Show All Patients

Patient Name Search

Exact Name (first last)

Patient ID Search

Image Comment Search

Search by Tooth Number

Date Range Search

Start Date

End Date

Dentist Name Search

Default Dentist (please update)



Double click a patient's name to see **each chart** that the patient has, organized by date. When you click the **specific chart date**, you will see all of the images in the **film strip** at the bottom

**Quick Search** lets you search by **preset search parameters**, like Yesterday/Today Studies

**Sub-Tabs** allow you to click on the **Patient Information** to see a patient's **entire information** entry

**Advanced Search** lets you search by **your own parameters**, like specific date ranges or by tooth #

### Patient Tab Notes:

1. Hover the mouse over any button in EVAsoft for the "tool tip."
2. Double click on a chart date and the chart will open up in the Chart View Tab.
3. Double click on any image in the film strip and the image will open in the Diagnosis View Tab.
4. The active patient is always named under the Patient Tab.
5. Delete studies for a patient in the Patient Information sub-tab.