

# eva ProIntegrate

## ProIntegrate Software Installation Instructions

### IMPORTANT NOTE

ProIntegrate software installation should be completed by a person with Windows O/S networking experience or by a trained professional.

### FILES

The following files make up the ProIntegrate product:

File Name	Description	Special Considerations
ProIntegrate.EXE	This is the main application executable.	Required file, must be installed.
ApteryxINI.DLL	This DLL contains the functionality required to handle the trial versions and registration of the product.	Required file, must be installed.
*.DRL	These various Data Drills contain the functionality required to obtain information from their respective database files. Examples: SoftDent.DRL; Dentrix.DRL	These Data Drills are the actual components that extract a patient's information from a 3 <sup>rd</sup> party database.
ProImage.DGDLL	This DGDLL file encapsulates the functionality that interfaces with the various ProImage installed objects used to open patients and exams in the ProImage application.	Required file, must be installed.
EVA.DGDLL	This DGDLL file encapsulates the functionality that interfaces with the various EVAsoft installed objects used to open patients and exams in the EVAsoft application.	Required file, must be installed.
ProIntegrate.CHM	This is the compiled on-line help for the ProIntegrate product.	This file is not mandatory, but without it, the ProIntegrate application will be unable to display its on-line help when the user clicks on the <b>Help</b> menu.

## HOW AND WHERE SHOULD I INSTALL PROINTEGRATE?

ProIntegrate should be installed on a central server, first, and then run (not installed) from each computer on the network. Install on the server before you run on the satellite computers. When you run ProIntegrate on the satellite computers (workstations) you must run ProIntegrate from the Server location. Where a network is not applicable, ProIntegrate can be installed on each computer in an office however this will require the purchase of a separate license for EACH computer.

## START MENU CONSIDERATIONS

When ProIntegrate runs it checks to see if it already exists in the user or the global Startup menu. If the application does not exist in either Startup menu, the user is prompted as to whether or not it should be added to the Startup menu. Dent-X recommends placing the shortcut to the ProIntegrate.EXE application in the global Startup menu.

## PROINTEGRATE SPECIFIC INFORMATION: WHAT IF THERE IS MORE THAN ONE DATABASE FROM THE SAME THIRD PARTY APPLICATION THEY WANT TO INTEGRATE WITH?

On occasion you may run into a situation where a user's site has more than one database from the same application that they want to bridge from (for example: a dual dentist office that runs two separate ACE Dental databases).

In the rare event that this occurs, the easiest way to handle this is to make a copy of the corresponding Data Drill file (for example: have a Ace.DRL and Ace\_2.DRL file in the application directory). The next time the application is started there will be two ACE Dental entries in the ProIntegrate options. Each entry can be set up to point to a different database.

---

## INSTALLATION NOTE

When installing the application it is recommended that you permit the user to select which applications they have installed instead of simply installing all of the ProIntegrate Data Drills at the same time. An example of how this is accomplished is included in this section.


Please make sure to install on the server first, before you run the ProIntegrate program from the workstations. Once ProIntegrate is installed successfully on the server run (open) the application, then continue to run (open) ProIntegrate on the workstations from the server location. At every workstation, when you run the software, you must run the ProIntegrate software from the server location.

For ProIntegrate to function the workstation must have the Practice Software and EVAsoft or ProImage already installed.




## INSTALLATION


1. Load the ProIntegrate CD into the CD-ROM drive of the server-computer. You must install on the server before you install anywhere else. The software should auto-run. Click "Next" when the following window appears:




2. Click "**Yes**" on the License Agreement window:




3. Click **“Next”** on the Destination Location window:




4. Select the Practice Management Software that will integrate with EVAsoft or ProImage:



5. Confirm that the correct Practice Software was selected then click **"Next"**




6. ProIntegrate will then install and it will let you know when the installation is complete. Click **"Finish"**





On the CD. Double click on the file DataGrabberProIntgrateBridgeUpdate.exe. Follow installation instructions. Additional upgrades may be available in the upgrade folder on this CD or at <http://www.apteryxware.com/datagrabber/downloads> .

7. Launch ProIntegrate by clicking in the **“Start”** menu, then **“All Programs”** menu, then **“ProIntegrate”** file, then the **“ProIntegrate”** application:



ProIntegrate will launch for the first time. Follow the prompts to enter in the registration information. An EVA icon will appear in the PC System Tray, next to the clock in the lower right-hand corner of the monitor. For Server 2008 & Windows 7 systems the icon may be hidden. Click on the up arrows in the system tray near the clock to show hidden icons and select customize. Then from the dialog select DataGrabber Application ProIntegrate (EVA Icon) and select Show Icon and notifications from the drop down menu.





## REGISTRATION

When you first run the application you are required to register ProIntegrate. Your registration number is contained on the ProIntegrate Product Registration Card in the CD case. Enter your registration number, as well as all other pertinent office information, into the appropriate text boxes on the "Register" pop-up window. If the remote server cannot be connected to, go to [www.apteryx.com/register](http://www.apteryx.com/register) to complete the registration process.

An example of the ProIntegrate Product Registration Card is:

<b>eva ProIntegrate</b>		
Practice Management Software Link for EVA and ProImage		
<table border="1"><tr><td>EVA ProIntegrate Registration Key</td></tr><tr><td>0123-4567-8901-2345</td></tr></table>	EVA ProIntegrate Registration Key	0123-4567-8901-2345
EVA ProIntegrate Registration Key		
0123-4567-8901-2345		
Registration is required for this product after the first 30-days If you encounter a problem registering please call 914-592-6665		

If you do not have a valid registration number you must click on **Continue Using Trial** button to activate the application as a 30-day trial and call Dent-X Technical Support.

## REGISTER USER

Once you enter in the ProIntegrate Product Registration number in the **Registration Number** field, enter in the **User Name**, **User Company**, **User Email** and **User Phone**

**Register**

# eva ProIntegrate

## STEP 1: REGISTRATION INFORMATION

Please enter the following information in order to register this product. If you were not already supplied with a registration number, please contact your appropriate technical support department at AFP Imaging.

Registration Number:

Install ID:

User Name:

User Company:

User Email:

User Phone:

## STEP 2: ACTIVATION CODE

Get Activation Code

You may obtain your Activation Code by clicking on the 'Get Activation Code' button if you have internet access. If you do not have internet access, please contact your appropriate technical support department at AFP Imaging.

## STEP 3: CONTINUE

Once all of this information is entered you must click the "Get Activation Code" button if the computer is online. If the computer cannot get online but the office has internet access on another computer then record the registration number and install ID and goto <http://reg.apteryx.com/cgi-bin/reg.pl> on the internet connected computer and register. After registration carefully record the activation code and type it exactly into the activation code box in ProIntegrate's registration dialog.

\*\*\* Once the Activation Code is entered, ProIntegrate is fully installed on the server. You will notice a small EVA ProIntegrate icon that will load into the PC System Tray beside the clock in the lower right hand of the computer monitor:






**IMPORTANT:** The ProIntegrate installation location on the server must be "shared" to access for networked workstations. If you do not know how to do this consult your local network administrator or IT support.

## INSTALLING ON A WORKSTATION


NOTE: For every workstation (satellite computer) do not install the software. Run the software, through the workstation, from the server location. To do this:

1. Double click "**My Computer**" to open it, then open the "**Network Drive**" (this is a different drive letter for every computer).
2. Open the folder "**Program Files**" or "**Program Files (x86)**" for **64 bit OS** by double clicking it.
3. Open the folder "**AFP Imaging**" by double clicking it.
4. Open the folder titled "**ProIntegrate**" by double clicking it.
5. Double click on the "**ProIntegrate**" file, which has an EVA icon next to it. This is the software that is located on the server. The icon will load into the system tray next to the clock. For Server 2008 & Windows 7 systems the icon may be hidden. Click on the up arrows in the system tray near the clock to show hidden icons and select customize. Then from the dialog select DataGrabber Application ProIntegrate (EVA Icon) and select Show Icon and notifications from the drop down menu. Optionally it is recommended that you create a desktop shortcut to the prointgrate.exe program.




6. Close all screens and *right click* on the  icon in the system tray.
7. Hover over the menu item "**Customize...**" then left click on "**Launched Applications....**" when it appears. Click continue.

8. Select either EVASoft if you are running EVASoft or ProImage if you are running ProImage and click OK.




9. Right click on the **eva** icon in the system tray.


10. Hover over the menu item "**Customize...**" then left click on "**Data Drills...**" when it appears.




11. If you get a warning here, click "**Continue**". This is fine.
12. Next, **double click on the data drill** listed, which is in accordance with what type of Practice Management software EVASoft or ProImage will link with.



13. Next, specify the path. This means that you will tell ProIntegrate where the Practice Management software is located (usually on the server also) so ProIntegrate will know where to look. Click on the button to the right in order to browse.




14. Next, "**Look in**" the **Network Drive** (whatever the corresponding drive letter is for that office) and double click on the folder for the Practice software.




15. Double click on the Practice software's **Data Folder**, and look for a file named **pat\_dat.dat**. Double click on this file and ProIntegrate will be mapped.

16. Repeat this procedure on every workstation, then you will be finished.

FINISH

Now ProIntegrate and  the icon is fully loaded and networked.

## USING PROINTEGRATE


ProIntegrate is to be used with either EVAsoft or ProImage open (that is running), as a software integration tool. While a user has a patient record selected in the Practice Management Software, simply click the EVA  ProIntegrate System Tray icon:

(Note: EVAsoft or ProImage should be added to the startup folder to ensure that it is opened and running .)


A series of prompts will appear notifying the user that the patient does not currently exist in EVAsoft or ProImage's database and will be added. The patient will be added to the ProImage database with all of the current information, according to the patient's record in the Practice Management Software.


## ADDITIONAL SERVICES

Try **right clicking** the EVA ProIntegrate icon to see a list of additional options and services:



Note: When ProIntegrate launches a pop-up window will appear:



Select **"Start Automatically When Windows Starts"** and click **OK**. Now the EVA ProIntegrate icon will always be available in the PC System Tray next to the clock in the bottom right hand  corner.



If the EVA ProIntegrate icon is ever not in the System Tray, then the user can simply click on the **“Start”** menu, then **“Programs”** menu, then **“ProIntegrate”** file, then the **“ProIntegrate”** application. The icon should load into the System Tray at this time. If the icon still cannot be seen then re-boot the computer to restart the program.

