



ImageWorks

Generations of Imaging

Computer* requirements for ImageWorks Panoura 18S with EVAsoft

ImageWorks strongly recommends that you purchase a new computer* for use with your Panoura 18S Digital Panoramic X-ray Scanner.

NOTE: This computer has several not typical requirements, double check the details of any computer you order.

- * Computer is not provided.
- ** Specifications are subject to change without notice.

- **Required Computer Specifications for ImageWorks Panoura with EVAsoft:**
 - Operation System - Microsoft Windows 7 **32 bit Professional**
 - Intel Quad Core 3.0 GHz
 - RAM - 4 GB (DDR3)
 - Available **PCI Express x1 slot.**
 - For dedicated Ethernet connection to the Panoura. This is separate and additional connection to any Local Area Network (LAN) connections.
 - A **Serial RS-232 Port¹** (this is the 'old style' 9 pin D sub connector, USB to RS232 adapters will **NOT** work).
 - CD-R / DVD-R/RW
 - Hard Drive - 250 GB
 - *Monitor with 1000:1 Contrast Ratio and 0.255 dot pitch*
 - Minimum resolution of 1280 x1024 @ 32 bits color depth.

Notice:

- ❖ Internet access is required for remote Technical Support assistance.
- ❖ It is strongly recommended that you consult your IT professional for a Data Back-up system.

¹ Many Desktop and most Notebook computer models do not include an RS232 port. If so an add-on on card will be required .