

Remote keyboard for Endos DC

The regular keyboard of Endos DC is a disc directly mounted on the timer cover.

The system can easily be transformed into a remote configuration, simply by adding a specific adaptation kit.

Using the adaptation kit, the keyboard which is originally installed on the timer cover can be detached and installed in a remote location by means of mounting brackets and accessories.

The "hole" in the timer cover can be covered with a blind adhesive sheet having the same design of the original keyboard.

The keyboard is connected to the main unit by means of a 15m (50ft) cable included in the adaptation kit.

The adaptation kit includes:

- Blind adhesive cover for the timer
- Mounting brackets for the remote keyboard
- 15m (50ft) Connection cable
- mounting accessories

This solution allows you to keep in stock some adaptation kits to be used whenever you need a remote system, without the need to keep different systems.

Standard mount

The keyboard is mounted directly on the timer box



Remote mount

The keyboard is detached from the timer box and installed outside the room



The Timer box is closed with a blind adhesive cover