



Prolmage

Prolmage Dental Image Management software meets the imaging demands of today's digital practices. Easy to learn and simple to operate, Prolmage has all of the advanced features necessary to deliver rapid, effective diagnoses. From the simplicity of navigation to the outstanding value it offers, Prolmage Dental Image Management software is the choice for any practice!

The Prolmage Experience:

- ✓ **Easy to Use - Simple Icon Navigation**
- ✓ **Automatic X-Ray Sequence Selection**
- ✓ **Automatic Image Save Function**
- ✓ **Endodontic Exposure Mode**
- ✓ **Decay and Abscess Detection Tool**
- ✓ **Manual and Automatic Image Contrast Enhancement**
- ✓ **Multiple Tooth Charts**
- ✓ **Send X-Rays to Insurance**

Prolmage Software:

- ✓ **EVA Digital X-Ray Sensor, Digital Pan & Intraoral Camera Software Suite**
- ✓ **Unlimited Multi-User Software per Network**
- ✓ **Stand-Alone or Network Environments**
- ✓ **Integrates with All Major Practice Software using ProIntegrate**
- ✓ **Developed and Supported in America**

The Value of Prolmage:

- ✓ **Free Lifetime Technical Support**
- ✓ **Free Workstation Licenses for Multiple Computer Networks**
- ✓ **Supports Almost All Intraoral Cameras**
- ✓ **Free Camera Module Included**
- ✓ **Free Cosmetic Imaging Module Included**

Recommended Computer Specs

Operating System: Microsoft 7 (32-bit), Vista Ultimate or Business (32-bit), Windows XP Professional (SP3 32-bit). Processor: Core 2 Duo (equivalent or better), RAM - 2 GB. General: USB 2.0 Ports, DVD-RW, 150 GB Hard Drive, High Quality Display (Flat Panel or CRT) with high brightness and DVI connection, LCD Monitor with 1000:1 Contrast Ratio and 0.255 dot pitch, Minimum resolution of 1600x1200 @ 32 bits color depth.

