

# MEDI-PRO/MP™ X-RAY FILM PROCESSOR

- Microprocessor Controlled
- Heavy Duty Rugged Construction
- Up to 200 Films Per Hour
- 90 Second Leading Edge Transport Time
- Exclusive Roller Transport System
- Easy Operation
- Attractive Stylish Design



VISIONARY IMAGING



Visionary Imaging

250 Clearbrook Rd.  
Elmsford, NY 10523  
Phone: 800.592.6666  
914.592.6100  
Fax: 914.592.6769  
e-mail: marketing@afpimaging.com  
www.afpimaging.com

AFP ... reliable and affordable medical imaging products for over two decades.

An ISO9001 Company

Automatic Film Processors

Veterinary Dental X-Rays

Digital Imaging

Panoramic X-Rays



Medi-Pro/MP control panel with developer and dryer temperature displays

Dealer:

**B**uilt for performance, the microprocessor controlled **Medi-Pro/MP™** X-ray film processor provides crisp, beautifully finished radiographs in just 90 seconds. A modular, compact design makes the Medi-Pro/MP an integral part of any facility's asset management program.

The Medi-Pro/MP is designed with great emphasis on structural integrity. With stainless steel construction and a high impact composite exterior, the Medi-Pro/MP is both robust and attractive. Precision machined, friction-reducing roller bearings, matched with an oversized drive motor enable the Medi-Pro/MP to transport film effortlessly. This results in ease of operation, low maintenance, and extended product longevity.

### Medi-Pro/MP Processor Specifications:\*

- Processing Time:** 90 seconds (leading edge transport time)
- Processing Capacity:** 51 in (129 cm) per minute  
200 14" X 17" (35cm X 43 cm) films per hour  
350 mixed film sizes per hour
- Dimensions:**
  - Processor:** Height: 21-1/8 " (53.6 cm)  
Width: 23" (58.4 cm)  
Length: 40-5/8" (103.2 cm)
  - Base:** Height: 24" (61 cm)
- Water Consumption:** 1/2 gal (1.9 liters) per minute - in processing mode  
(No water consumption in standby mode)
- Tank Capacities:** Developer: 2-3/4 gal (10.4 liters)  
Fixer: 1-3/4 gal (6.6 liters)  
Wash: 1-3/4 gal (6.6 liters)
- Weight:** 250 lbs (114 kg) - without fluids
- Power Requirements:** 120V, 20A, 50/60 Hz  
240V, 10A, 50/60 Hz
- Drain Requirements:** 4 gal (15 liters) per minute (required only when tanks are emptied during periodic maintenance)
- Warranty:** 2 year limited on parts
- Part Numbers:** Medi-Pro/MP processor 9992312100 (120V)  
Medi-Pro/MP processor 9992312200 (240V)  
**Note:** A base stand and two (2) 7-gal replenishment tanks are included with both part numbers listed above

\*Specifications are subject to change without notice

