

NewTom VGi Flex

Mobile Cone Beam 3D Imaging System



NewTom and ADI

Own Your Own Mobile Clinic

NewTom VGi FLEX can be purchased in a variety of configurations. From a fully outfitted, turn-key mobile imaging operation, to a NewTom VGi Flex by itself, we offer a number of options for setting up your mobile imaging business. AFP teamed up with mobile solutions provider ADI Mobile Health to supply a “Standard” option, and a more “Well-Equipped” version of a mobile dental radiology center. Both versions, built on a Dodge Sprinter chassis, are ready to begin scanning patients with just the turn of the key. Our expert installation team can also provide the requirements necessary for placing NewTom VGi Flex in the van or truck of your choice

Standard VGi Sprinter \$85k*

- 2009 Dodge Sprinter 3500 Chassis 4x2 170”
- Wide Body/Flat Floor – LHD, Interior Dimensions: L-170”,W-94”, H-85”,
- V6 Turbo Diesel Engine, 5 Speed Automatic, with Dual Airbags
- AM/FM/CD/MP3
- Remote door locks, Stability control, Power windows

*does not include the cost of NewTom VGi Flex



- ABS, 4-wheel disc brakes
- Refrigerated/cooled storage box
- 4 way adjustable seats
- Exhaust & Heat Riser System
- Ducted AC, w/heat strip and thermostat
- Recessed Fluorescent lighting
- Fiberglass Insulated interior/ exterior walls and acoustic ceiling
- Auxiliary battery compartment, isolated from engine battery
- Right rear entrance door with window and manual step
- Basement storage compartment
- Exterior Porch light
- Computer Desk and Laminated cabinets
- 3D Mounting wall
- Heavy Duty Vinyl Flooring
- 5 120V interior outlets
- 120V AC service Panel
- 12V DC panel with breakers
- 12V DC charger/converter
- 3 Fire Extinguishers

Well-Equipped VGi Sprinter \$95k* (Standard VGi Sprinter plus the following):


- Shaded windows
- Back-up & Side view Camera System
- Back-up alarm
- Side to side Stabilizing/Leveling Jacks
- Air Suspension
- Additional AC/ heater
- 3 passenger bench seat in cab
- Broadband cellular fax modem with Wi/Fi

*does not include the cost of NewTom VGi Flex


Other Options

- 8 foot Power awning
- Landline power cord with Surge Guard
- Carpeted floor instead of vinyl
- Overhead cabinets with lights (Desk or Bench)
- Wheel Chair Ramp
- Automatic steps
- Bi-Fold X-ray Door
- TV with DVD player
- 1000 Watt Inverter (UPS) with 2nd Battery
- Toe kick or wall electric heaters

Sample Standard Plan



Sample Well-Equipped Plan



Manufactured in Italy by QR srl,
an AFP Imaging Company
www.afpimaging.com

AFP Imaging Corporation
250 Clearbrook Road
Elmsford, NY 10523 USA

1 800.592.6666
914.592.6100