

Put 3D Diagnostic Imaging Capability in Your Office Today!



Available immediately, Certified Pre-Owned
NewTom 3G's with 12" FOV

AFP Imaging has a limited available inventory of certified, pre-owned NewTom 3G CB3D scanners. These units were installed at a reputable, highly-successful dental imaging center, and have been maintained and recently upgraded by factory trained engineers.

All machines are capable of 6", 9", and 12" FOV's and have been used by the imaging center for everything from implants to orthodontics and maxillofacial surgery cases, as well as Sleep Apnea and Airway applications.

With NewTom 3G you get:

- Improved 3D from upgraded cameras
- Latest Dual Core Computers = less than 1 min of processing time
- Voxel Sizes:

Voxel Size Matrix	12"	9"	6"
Small Field	.38mm	.25mm	.18mm
Large Field	.42mm	.29mm	.21mm

- Compatible with 3rd Party Software packages



- 1 Year warranty included (Extended Service Agreements Available)
- Includes 2 Day training by Imaging Center employee
- 40 hours of complimentary technical support from the Imaging center


See examples of the image quality you can expect from these machines at www.cdental.com.

For more information, please contact AFP Imaging (800) 592-6666

OR

Ext 101 for Daniel Fields—dfields@afpimaging.com

Ext 217 for Michael Ellison—mellison@afpimaging.com



AFP Imaging
250 Clearbrook Road
Elmsford, NY 10523
(914) 592-6100
www.afpimaging.com

