

# ImageV™ PACS

The ImageV™ PACS is a digitally integrated turn-key solution developed specifically for the veterinary market that incorporates full body and dental acquisition in one central workstation. This powerful system allows customers utilizing AFP veterinary products such as EVA® Vet and the DIGIVET® digital product series, to integrate for acquiring, manipulating and storing images. Veterinary practices now need only one software program to perform these tasks. ImageV PACS makes digital imaging and diagnostic viewing easier than ever before.

With ImageV PACS, veterinary viewing for CR, CCD, DR, ultrasound, endoscopy, and digital dental modalities is made simple. The system's sophisticated DICOM-compliant PACS server can accept virtually any DICOM image and export it in a variety of formats to additional viewing stations, multiple exam rooms, CD or DVD, DICOM send, or remote locations.

AFP's market-leading EVA Vet can capture directly into the ImageV Dental PACS module as well. Within the Dental Module, you can choose or create a dental mount template and acquire images at the click of a button. Dental and Whole Body images are stored together, under one patient.

ImageV PACS can store anywhere from 10,000-100,000 digital images and is developed to grow with the clinic and adapt with your practice as needed. The ImageV PACS is the only veterinary PACS system offering a robust archiving system and open architecture.



## Why ImageV™ PACS


- Images are instantly sent to a PACS archive server, then distributed to viewing stations
- Integrate seamlessly with Practice Management Software
- Maximize workflow with seamless integration for capturing, archiving and distributing images throughout the practice
- Burns CDs and/or DVDs for owners, clients and referring doctors
- Provides data security
- Reduced wait/response time
- Enhanced patient care
- One computer serves all digital imaging needs
- Review images quickly and accurately with the user-friendly interface

# ImageV™ PACS Advantages

Can utilize all AFP Digital offerings and the ImageV PACS technology all in one system	Dental software can be integrated for a complete full body/dental solution
Incorporates full body and dental acquisition into one DICOM work list	Switch from one modality to the next with one-touch Plug and Play functionality
Easily manage and store your medical records with high volume storage capability	Seamless Practice Management software (PMS) Integration with ImproMed®, IntraVet and AVImark®
One vendor installing digital system, avoids confusion and communication delays	Critical for disaster recovery
Expands with your practice	Data remains proprietary to customer
One workstation does it all!	

SMO-157-00

Copyright © 2009 AFP Imaging Corporation. All rights reserved. Reproduction of this whole or any part of the contents without written permission is prohibited.  
EVA Vet, ViewAll and DIGIVET: Registered Trademarks of AFP Imaging Corp. Printed in the U.S.A.



## Call us for expert digital solutions and support.

The Total Veterinary Imaging Solution from AFP Imaging

EVA® Vet Dental Digital Sensors | Image-Vet 70ACP Dental X-Ray | DIGIVET® CR and CADR  
Digital VetTek CCD Tables | DIGIVET® Equine DR | ViewAll® Vet | ImageV™ PACS



**AFP Imaging Corporation**  
250 Clearbrook Road  
Elmsford, NY 10523 USA

1.800.592.6666  
P 914.592.6100  
F 914.592.6148

[www.afpimaging.com](http://www.afpimaging.com)

**Total Imaging Solutions | Veterinary | Dental | Medical | NewTom Cone Beam 3D**